



Disposable Linear Cutter  
for Endoscope Use

## 一次性使用旋转切割吻合器

注册号：浙食药监械（准）字2014第2080610号

浙江格宜医疗器械有限公司

ZHEJIANG GEYI MEDICAL INSTRUMENT CO.,LTD.

## 一次性使用旋转切割吻合器

Disposable Linear Cutter for Endoscope Use

### 产品特性:

- 1.组件可左右45度旋转, 灵活到达不同手术部位
- 2.三排钉高不同, 确保吻合效果
- 3.扳机与手柄之间的小跨度, 增加了击发时的舒适度, 减小击发用力
- 4.击发时的触觉反馈, 可使外科医生确认缝钉已成型
- 5.器械杆可以360度旋转, 使器械可以在任意部位使用
- 6.直观、简单的手柄设计, 可使医生一目了然, 操作方便

### 适用范围:

适用于腔镜和开放手术的消化道组织离断、切除和吻合。

### Application Scope:

Suitable for dividing, resection and anastomosis in the endoscopic surgery or open surgery.

### Product characteristics:

- 1.Components can be rotated 45 degree left and right, flexible to reach different surgical site.
- 2.Three rows of staple height are different, which ensure the consistent results.
- 3.The small span between the trigger and handle increases the fires of comfort and reduces the firing strength
- 4.The tactile feedback when firing, allows the surgeon to confirm that suture staple has shaped
- 5.Instrument can be rotated 360 degree, which makes the equipment can be used in any place.
- 6.Intuitive and simple handle design makes the doctor be clear at a glance, easy to operate.

### 规格型号 Model & Specification :

旋转切割吻合器型号 Stapler Model	切割组件型号 Cutting Component Model	缝合长度(mm) Suture Length(mm)	杆部长度(mm) Rod Length(mm)	钉高 Staple Height	缝钉数(个) Staple Quantity(pcs)
GYRC-60	GYRCZ-(A/B)30(V/L/H/M)	31	60	V:2.0	42
	GYRCZ-(A/B)45(V/L/H/M)	46		L:2.5	66
	GYRCZ-(A/B)60(V/L/H/M)	61		H:3.5 M:4.8	90
GYRC-100	GYRCZ-(A/B)30(V/L/H/M)	31	100	V:2.0	42
	GYRCZ-(A/B)45(V/L/H/M)	46		L:2.5	66
	GYRCZ-(A/B)60(V/L/H/M)	61		H:3.5 M:4.8	90
GYRC-160	GYRCZ-(A/B)30(V/L/H/M)	31	160	V:2.0	42
	GYRCZ-(A/B)45(V/L/H/M)	46		L:2.5	66
	GYRCZ-(A/B)60(V/L/H/M)	61		H:3.5 M:4.8	90
GYRC-260	GYRCZ-(A/B)30(V/L/H/M)	31	260	V:2.0	42
	GYRCZ-(A/B)45(V/L/H/M)	46		L:2.5	66
	GYRCZ-(A/B)60(V/L/H/M)	61		H:3.5 M:4.8	90